



2010  
Good Design Award



Oxygen concentrator is a device that generates oxygen with a purity of  $92 \pm 3$  by using the environment air that comprehends in the ratio of %21. Oxygen purification system designed appropriate to the human health. Today, many people are living with different illnesses, such as, heart, lung and the other respiration diseases. Most of these people have found help from the oxygen treatments that they take at home or hospital.

This device will help you to take adequate oxygen treatment that you need.

Using the oxygen without permission of authorized people can be dangerous.

Device should be used under the supervision of a licensed physician after reading the instructions manual.

### SPECIFICATION

OXYGEN FLOW RANGE	1- 5L / min
OXYGEN CONCENTRATION	%92 $\pm$ 3
MAXIMUM OUTLET PRESSURE	8 PSI $\pm$ 1
VOLTAGE	230 VAC $\pm$ %10, 50 Hz
CURRENT	1,40 AMP $\pm$ %10
OPERATING TEMPERATURE RANGE	10°C - 40°C
OPERATING HUMIDITY RANGE	10 - 90 % RH
WEIGHT	14 kg
NOISE LEVEL	< 45 dB

### ACCESSORIES

