



FLAT REELS

Flat reels are composed of one side medical kraft paper, one side medical plastic film. The raw materials used are certified and offer a first-class seal, free from paper residues or film lamination and a constant welding during the sterilization.

GUSSETTED REELS

Gussetted reels are composed of one side medical kraft paper, one side medical plastic film. The raw materials used are certified and offer a first-class seal, free from paper residues or film lamination and a constant welding during the sterilization. The gussetted allows the packing of cumbersome products.




FLAT POUCHES

Flat pouches are composed of one side medical kraft paper, one side medical plastic film. These pouches are suitable for Ethylene Oxide and Steam Sterilization.

ECS TAPES

The adhesives ECS tapes are available in many sizes and heights with or without the indicator, suitable for steam, dry heat and ETO sterilization.



CREPED PAPER

High quality creped paper for sterilization packaging, TNT NON-WOVEN for sterilization packaging and BLUE SMS.